



Quick Reference Environmental Cleaning for School Nurses

How to clean equipment and surfaces, when to clean, what to use

Visibly soiled surfaces must be cleaned with friction first then re-wiped. Non-visibly soiled surfaces may be wiped thoroughly with friction once. When using spray disinfectant, spray/saturate towel and wipe surface with friction. Spraying the surface and letting it dry is not adequate to clean and disinfect. Wear gloves.

| Item to clean | Wipe down what? | When? | With what? |
|---|--|--|---|
| Blood pressure cuff | Cuff, tubing, bulb (if manual) | After each use | Wipe with low-level disinfectant (LLD); if visibly soiled, wash in soap and water, rinse and hang to dry. |
| Stethoscope | Bell and tubing | After each use | Wipe with IPA (isopropyl alcohol) or LLD. |
| Reflex hammer | Handle and head | After each use | Wipe with IPA or LLD. |
| Otoscope speculae (Disposable preferred) | If reusable, wash and disinfect speculae | All surfaces after each use | IPA may be used for non-disposable otoscope speculae, soak for 20 minutes. |
| Otoscope handle | Handle | All surfaces after each use | Wipe with LLD and air dry. |
| Ear thermometer, Temporal scanner thermometer | All surfaces | All surfaces after each use | Wipe with LLD and air dry. |
| Digital thermometer with disposable probes | All surfaces, except disposable probe | After each use | Wipe "machine" with LLD and air dry; Discard disposable probe after each use. |
| Basin (Disposable preferred) | Basin | After each use | Wash in enzymatic detergent and rinse well, then autoclave. <i>(If do not have autoclave, use disposable-one time use only)</i> |
| Canes, walkers, crutches, wheelchairs, rehab equipment | Special attention to surfaces that are in contact with people | Between patients | If visibly soiled, clean first with friction. Then wipe down with LLD. |
| Exam tables, gurneys | Top of table, other areas that came in contact with patient and/or bodily fluids | After each use | If visibly soiled, clean first with friction. Then wipe down with LLD. |
| Door knobs, phones, keyboards, light switches, hand sanitizer dispenser, pen lights & other "hand touch" items in nursing area | "Community" pen at the desk; Shared keyboards, counters, telephones, doorknobs, drawer pulls and other "hand touch" areas | At least twice daily <i>(Clean often because this is a "healthcare" area)</i> | If visibly soiled, clean first with friction. Then wipe down with LLD. Disinfect keyboards for 5 seconds daily and when visibly dirty by wiping with LLD or IPA. |
| Waiting room: chairs, tables, etc. | All surfaces that can come in contact with patient | First thing in the morning and/or at the end of the day | If visibly soiled, clean first with friction. Then wipe down with LLD. |

Guide to Sterilization and Disinfection of Equipment & Surfaces

You MUST thoroughly clean items to remove all visible soil before disinfecting or sterilizing

| Intended Use | Level of Process Required | Products* *Brand names are used as examples only, no endorsement is implied. | | | | | | | | | | | | | | |
|--|--|--|------------------|----------------------|----------------------------|--------|-------------------------------------|----------------------------|------------------------|------|----------------------------|------------|----------------------------|------------|-------------------------------|------------|
| <p>Critical</p> <p>Intended use: Objects that enter normally sterile tissue, the vascular system or through which blood flows, should be sterile (<i>instruments</i>)</p> | <p>Sterilization</p> <p>(Steam, gas, hydrogen peroxide plasma, or chemical sterilization)</p> | <p><u>For chemical sterilization:</u></p> <p>Glutaraldehyde (≥ 2.0%) (Cidex, Metricide) Hydrogen peroxide – HP (7.5%) (Sporox) Peracetic acid – PA (0.2%) HP (1.0%) & PA (0.08%) HP (7.5%) & PA (0.23%) Glut (1.12%) & Phenol/phenate (1.93%) <i>(Exposure time on all per manufacturers' recommendations)</i> HP is NOT your clinic wound cleaning product!</p> | | | | | | | | | | | | | | |
| <p>Semi-critical</p> <p>Objects that touch mucous membranes or skin that is not intact require an HLD process (<i>scissors</i>)</p> | <p>High Level Disinfection (HLD)</p> <p><i>(FDA regulates products)</i></p> | <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; border-bottom: 1px solid black;"><u>Germicide</u></th> <th style="text-align: right; border-bottom: 1px solid black;"><u>Concentration</u></th> </tr> </thead> <tbody> <tr> <td>Glutaraldehyde</td> <td style="text-align: right;">≥ 2.0%</td> </tr> <tr> <td>Ortho-phthalaldehyde (OPA) (12 min)</td> <td style="text-align: right;">0.55%</td> </tr> <tr> <td>Hydrogen peroxide (HP)</td> <td style="text-align: right;">7.5%</td> </tr> <tr> <td>HP and paracetic acid (PA)</td> <td style="text-align: right;">1.0%/0.08%</td> </tr> <tr> <td>HP and PA</td> <td style="text-align: right;">7.5%/0.23%</td> </tr> <tr> <td>Hypochlorite (free chlorine)*</td> <td style="text-align: right;">650-675ppm</td> </tr> </tbody> </table> <p>*May cause cosmetic and functional damage <i>(Exposure time ≥ 12 min to 30 min @ 20° C, see manufacturers' recommendations)</i></p> | <u>Germicide</u> | <u>Concentration</u> | Glutaraldehyde | ≥ 2.0% | Ortho-phthalaldehyde (OPA) (12 min) | 0.55% | Hydrogen peroxide (HP) | 7.5% | HP and paracetic acid (PA) | 1.0%/0.08% | HP and PA | 7.5%/0.23% | Hypochlorite (free chlorine)* | 650-675ppm |
| <u>Germicide</u> | <u>Concentration</u> | | | | | | | | | | | | | | | |
| Glutaraldehyde | ≥ 2.0% | | | | | | | | | | | | | | | |
| Ortho-phthalaldehyde (OPA) (12 min) | 0.55% | | | | | | | | | | | | | | | |
| Hydrogen peroxide (HP) | 7.5% | | | | | | | | | | | | | | | |
| HP and paracetic acid (PA) | 1.0%/0.08% | | | | | | | | | | | | | | | |
| HP and PA | 7.5%/0.23% | | | | | | | | | | | | | | | |
| Hypochlorite (free chlorine)* | 650-675ppm | | | | | | | | | | | | | | | |
| <p>Non-critical</p> <p>Objects that will not come in contact with mucous membranes or non-intact skin (e.g. environmental surfaces) require a low level process that kills vegetative bacteria, fungus and some viruses (Hepatitis B, C, MRSA and HIV).</p> | <p>Low Level Disinfectant (LLD)</p> <p><i>(EPA regulates hospital-level products)</i></p> | <table style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="text-align: left; border-bottom: 1px solid black;"><u>Germicide</u></th> <th style="text-align: right; border-bottom: 1px solid black;"><u>Concentration</u></th> </tr> </thead> <tbody> <tr> <td>Ethyl or isopropyl alcohol</td> <td style="text-align: right;">70-90%</td> </tr> <tr> <td>Chlorine</td> <td style="text-align: right;">100ppm (1:500 dilution) **</td> </tr> <tr> <td>Phenolic</td> <td style="text-align: right;">*</td> </tr> <tr> <td>Iodophor</td> <td style="text-align: right;">*</td> </tr> <tr> <td>Quaternary ammonium (quat)</td> <td style="text-align: right;">*</td> </tr> </tbody> </table> <p>*Use manufacturers' recommendations for concentrations **5.25-6.15% household bleach diluted 1/500 provides > 100ppm available chlorine <i>(Exposure time ≥ 1 min)</i></p> <p>EPA registered, pop-up towelettes are usually quats and are effective low-level disinfectants.</p> | <u>Germicide</u> | <u>Concentration</u> | Ethyl or isopropyl alcohol | 70-90% | Chlorine | 100ppm (1:500 dilution) ** | Phenolic | * | Iodophor | * | Quaternary ammonium (quat) | * | | |
| <u>Germicide</u> | <u>Concentration</u> | | | | | | | | | | | | | | | |
| Ethyl or isopropyl alcohol | 70-90% | | | | | | | | | | | | | | | |
| Chlorine | 100ppm (1:500 dilution) ** | | | | | | | | | | | | | | | |
| Phenolic | * | | | | | | | | | | | | | | | |
| Iodophor | * | | | | | | | | | | | | | | | |
| Quaternary ammonium (quat) | * | | | | | | | | | | | | | | | |

Note: Intermediate level disinfectants - usually used for therapy or whirlpool tubs.