

Division of Laboratory Services 630 Hart Lane Nashville, TN 37216 615-262-6300

https://www.tn.gov/health/health-program-areas/lab.html

Disease/Agent Suspected or Test Requested:

RPR (Rapid Plasma Reagin)—Syphilis

Provider Requirements	
Acceptable Specimen Sources/Type(s) for Submission	• Serum
TDH Requisition Form Number	PH-4182
Media Requirements	 Serum Red Top or Serum Separature Tube (SST) Sterile, plastic screw capped vial (serum)
Special Instructions	 PTBMIS Order Code: 86592 RPR (Syphilis Test) StarLIMS Order Code: 3302 (Syphilis Panel) Specimens may be stored for up to 7 days at 2-8°C. If testing will be delayed more than 7 days for specimens stored at 2-8°C, aliquot serum, and store frozen at -20°C or colder.
Shipping Instructions	 Specimens should be shipped to the Laboratory the same day they are collected. Ship specimens on cold packs. If specimen has been frozen, ship on dry ice.
Laboratory Section Performing Testing	Serology
Lab Location(s) Performing Test	Nashville; Knoxville

All infectious substance shipments must conform to U.S. Department of Transportation (DOT) Hazardous Materials Regulations (HMR 49 C.F.R Parts 171-180).