



Division of Laboratory Services
 630 Hart Lane
 Nashville, TN 37216
 615-262-6300

<https://www.tn.gov/health/health-program-areas/lab.html>

Disease/Agent Suspected or Test Requested:

Tuberculosis – Quantiferon Gold

Provider Requirements	
Acceptable Specimen Sources/Type(s) for Submission	Whole Blood
TDH Requisition Form Number	Quantiferon TB Gold Log Sheet
Media Requirements	Quantiferon TB Gold blood collection tubes
Special Instructions	PTBMIS Order Code: 86480 (TB Blood Test Gold In-Tube) StarLIMS Order Code: 5008 (TB Panel)
Shipping Instructions	<ul style="list-style-type: none"> • Keep sample set together. • Incubate for 16-24 hours before shipping. • Must be received within 72 hours after incubation.
Laboratory Section Performing Testing	Mycobacteriology
Lab Location(s) Performing Test	Nashville

All infectious substance shipments must conform to U.S. Department of Transportation (DOT) Hazardous Materials Regulations (HMR 49 C.F.R Parts 171-180).