

Division of Laboratory Services 630 Hart Lane Nashville, TN 37216 615-262-6300

https://www.tn.gov/health/health-program-areas/lab.html

West Nile Virus - Encephalitis/Human

Provider Requirements	
Acceptable Specimen Sources/Type(s) for Submission	CSF Paired Serum Samples:
TDH Requisition Form Number	PH-4182
Media Requirements	 Red-stoppered vacuum tube (whole blood) Sterile, plastic screw capped vial (serum/CSF)
Special Instructions	 PTBMIS Order Code: WSTNILE (West Nile Virus) StarLIMS Order Code: 4001 (West Nile Virus Panel) Specimen should be kept refrigerated. If testing will be delayed more than 48 hours post collection, serum should be aliquoted and stored at -20°C or colder.
Shipping Instructions	If specimen has been refrigerated, ship on cold packs. If specimen is frozen, ship on dry ice.
Laboratory Section Performing Testing	Knoxville Regional Laboratory
Lab Location(s) Performing Test	Knoxville

All infectious substance shipments must conform to U.S. Department of Transportation (DOT) Hazardous Materials Regulations (HMR 49 C.F.R Parts 171-180).