



Division of Laboratory Services
 630 Hart Lane
 Nashville, TN 37216
 615-262-6300

<https://www.tn.gov/health/health-program-areas/lab.html>

Disease/Agent Suspected or Test Requested:

**Chlamydia trachomatis/Neisseria gonorrhoeae (Gonorrhea) antigen detection
 by Nucleic Acid Amplification (NAAT)**

Provider Requirements	
Acceptable Specimen Sources/Type(s) for Submission	<ul style="list-style-type: none"> • Endocervical • Urethral (male) • Rectal • Throat • Vaginal • Urine
TDH Requisition Form Number	PH-4182
Media Requirements	<ul style="list-style-type: none"> • Aptima Urine Specimen Collection Kit for Male and Female Urine Specimens • Aptima Multitest Swab Specimen Collection Kit for Rectal, Throat and Vaginal Swab Specimens • Aptima Unisex Swab Specimen Collection Kit for Endocervical and Male Urethral Swab Specimens
Special Instructions	<p>PTBMIS Order Codes:</p> <ul style="list-style-type: none"> • GENPU (GEN PROBE, URINE) • GENPP (GEN PROBE, UROGENITAL PENILE SW) • GENPE (GEN PROBE ENDOCERVICAL) • GENPO (GEN PROBE THROAT) • GENPV (GEN PROBE VAGINAL) • GENPA (GEN PROBE ANAL) <p>StarLIMS Order Code: 2199 (CT/GC Panel)</p>
Shipping Instructions	Ship specimens at room temperature or on cold packs.
Laboratory Section Performing Testing	Serology
Lab Location(s) Performing Test	Nashville; Knoxville

All infectious substance shipments must conform to U.S. Department of Transportation (DOT) Hazardous Materials Regulations (HMR 49 C.F.R Parts 171-180).